



Snap-In Syringe Holder

Holds 20, 30 and 60 mL Syringes

The new and improved Snap-In Syringe Holder is 3D printed and engineered to the tightest tolerances to fit securely in the patented rotating carousel in the RS 3400, providing improved durability, security and ease of use — and it has a more flexible fit for multiple syringe volume sizes.

The RS 3400 provides the convenience, flexibility and efficiency needed for today's busy processing runs.

Volume Flexibility, Product Security and Ease of Use

Improved Durability and Increased Life Span

- Made of carefully researched and selective Nylon material, the new Snap-In Syringe Holder provides the benefits of more flexibility, durability and increased life span.

Improved Syringe and Holder Fit


- 20, 30 and 60 mL syringes fit within the syringe holder without sticking out and causing interference.
- Latch grooves have been added to the outer wall of the holder to more securely grip the syringe holder in place within the canister holder's spring retention feature.

Improved Syringe Security

- Retaining and lock-in prongs hold the syringes secure and in place, negating the need for any rubber bands to keep the syringes secure within the holder.

Ordering Information

- Snap-In Syringe Holder: P/N 1427100
- Optional and can be purchased separately.



The Snap-In Syringe Holder works with legacy versions of RS 3400



4907 Golden Parkway, Suite 400
Buford, GA 30518
678-765-7900
RadSource.com

